


# Nitinella Plus

SX-ELLA Stent Biliary



# Nitinella Plus



Nitinol braided self-expanding stent

Silicone covering

Kink-resistant delivery system

Endoscopic and percutaneous access



## Basic description

### SX-ELLA Stent Biliary - Nitinella Plus

Nitinella Plus stent is designed for the palliative treatment of biliary strictures caused by malignant neoplasms



#### Innovative features / Benefits

- ☑ Excellent flexibility / Smooth placement in tortuous anatomies
- ☑ Radial force of stent is designed to maintain patency and resist migration
- ☑ Optimized stent shortening / Facilitated stent positioning
- ☑ Atraumatic stent ends / Effective prevention of trauma to the bile duct wall
- ☑ Reconstrainability / Enhanced physician control in deployment and repositioning
- ☑ Radiopaque markers at both stent ends and in the middle / Designed to aid stent placement
- ☑ Durable silicone covering / Prevention of tumor in-growth
- ☑ Braided kink-resistant delivery system / Superior system trackability and stent deployment

#### Additional description



Delivery system is tipped with radiopaque olive

Recommended guide wire: 0.035" (0.89 mm) ultra-stiff with minimum length of 360 cm (Endoscopic access) or 180 cm (Transhepatic access)



## Available sizes

Endoscopic delivery system - 175 cm						Transhepatic delivery system - 60 cm		
REF No.						REF No.		
Stent diameter [mm]	Stent length [mm]	Delivery system diameter [F]	Uncovered stent	Partially covered stent	Fully covered stent	Uncovered stent	Partially covered stent	Fully covered stent
8	40	7 for uncovered / 9 for covered stent	019-01B-08-040	019-01C-08-040	019-01D-08-040	019-01B-08-040-T	019-01C-08-040-T	019-01D-08-040-T
	60		019-01B-08-060	019-01C-08-060	019-01D-08-060	019-01B-08-060-T	019-01C-08-060-T	019-01D-08-060-T
	80		019-01B-08-080	019-01C-08-080	019-01D-08-080	019-01B-08-080-T	019-01C-08-080-T	019-01D-08-080-T
	100		019-01B-08-100	019-01C-08-100	019-01D-08-100	019-01B-08-100-T	019-01C-08-100-T	019-01D-08-100-T
10	40	7 for uncovered / 9 for covered stent	019-01B-10-040	019-01C-10-040	019-01D-10-040	019-01B-10-040-T	019-01C-10-040-T	019-01D-10-040-T
	60		019-01B-10-060	019-01C-10-060	019-01D-10-060	019-01B-10-060-T	019-01C-10-060-T	019-01D-10-060-T
	80		019-01B-10-080	019-01C-10-080	019-01D-10-080	019-01B-10-080-T	019-01C-10-080-T	019-01D-10-080-T
	100		019-01B-10-100	019-01C-10-100	019-01D-10-100	019-01B-10-100-T	019-01C-10-100-T	019-01D-10-100-T



Distributor:

Manufacturer:

**ELLA-CS, s.r.o.**

Milady Horakove 504/45, Trebes,  
500 06 Hradec Kralove, Czech Republic

Tel.: +420 495 279 111 Fax: +420 495 265 655  
E-mail: info@ellacs.eu URL: www.ellacs.eu